



Stand Out from the Crowd

Typical RP Columns
Similar selectivity limiting your choices for challenging separations



Platinum Columns
Unique selectivity that succeeds when traditional reversed phase columns fail



Platinum™ HPLC Columns — Controlled Silica Exposure

The Platinum Column Advantage

Controlled silica exposure is the difference that makes Platinum Columns unique. The approach taken by most column manufacturers in making reversed-phase columns is to thoroughly cover the silica with bonded phase to minimize any interaction between polar analytes and the packing medium's silica backbone. Platinum media takes a different approach. Instead of maximizing the coverage of bonded phase to hide the silica, the exposure of the silica is controlled to provide a dual mode separation medium with both polar and non-polar sites exposed to your samples. This extends polar selectivity well beyond what other reversed-phase columns offer and gives separations other columns cannot.


Standard Platinum Columns vs Platinum EPS Columns

Platinum Columns come in two varieties offering different levels of silica exposure. Standard Platinum has a moderate silica exposure and is best used with neutral and moderately polar compounds. Platinum EPS (Extended Polar Selectivity) has a high level of silica exposure and is best used with compounds containing more than two polar functional groups.

Most Reversed Phase Columns Have the Same Selectivity

Most reversed phase columns have the same selectivity. If one reversed phase column does not separate the sample, chances are other reversed phase columns will also fail!

Basic Pharmaceuticals On Competitive Reversed Phase Columns




1. Procaine
2. Lidocaine
3. Tetracaine
4. Diphenhydramine

Column: C18, 5µm, 100 x 4.6mm
Mobile Phase: Methanol:0.025M K₂HPO₄, pH 7.0 (85:15)
Flow Rate: 2.0mL/min
Detector: UV at 230nm

The Platinum Column Difference

Platinum Columns produce results when other columns cannot. Platinum's controlled silica exposure extends the selectivity range well beyond that of other reversed phase columns.

Basic Pharmaceuticals



1. Procaine
2. Lidocaine
3. Tetracaine
4. Diphenhydramine


Column: C18, 5µm, 100 x 4.6mm
Mobile Phase: Methanol:0.025M K₂HPO₄, pH 7.0 (85:15)
Flow Rate: 2.0mL/min
Detector: UV at 230nm

The GRACE Davison Product Lines

Unique Reversed Phase Selectivity

Trying to solve difficult separation problems using typical reversed phase columns often leads to the same problem. Choose Platinum™ Columns for completely different selectivity. See chart below.

Hydrophobic/Hydrophilic Balance at pH 3 of Various Reversed Phase HPLC Packings



Plotting k' values of different compounds (polar vs. nonpolar) demonstrate the unique selectivity of Platinum & Platinum EPS Columns, compared to conventional reversed phase columns.

The GRACE Davison Product Lines

Alltech

DAVISIL

CROM

JONES
CHROMATOGRAPHY

MODcol

VYDAC

Polar Selectivity Benefits


Complementary Selectivity

Standard Platinum™ and Platinum EPS Columns produce different retention times and selectivities for the same sample.


Sedatives/Hypnotics Mix

CHROM
8187

Standard Platinum C18 Column



Platinum EPS C18 Column



NOTE
Standard Platinum Columns resolves an impurity peak that coelutes with Triazolam on Platinum EPS Columns.

1. Barbitol
2. Phenobarbital
3. Butalbital
4. Oxazepam
5. Triazolam
6. Unknown
7. p-Xylene

Columns: 5µm, 150 x 4.6mm
Mobile Phase: Methanol:Acetonitrile: 0.05M NH₄H₂PO₄ pH 4.5 (50:10:40)
Flow Rate: 1.0mL/min
Detector: UV at 230nm

CHROM
8188

Selectivity Values*

Compound	Standard Platinum C18	Compound	Platinum EPS C18
1. Barbitol	17.65	1. Barbitol	13.48
2. Phenobarbital	9.66	2. Phenobarbital	8.00
3. Butalbital	5.97	3. Butalbital	8.28
4. Oxazepam	2.80	4. Oxazepam	2.11
5. Triazolam	2.19	5. Triazolam	0.86
7. p-Xylene	k' ₁ =4.06	7. p-Xylene	k' = 2.48

*All selectivity values determined using p-xylene.

Symmetrical Peaks for Difficult Samples


Peak shapes for difficult amine- or carboxyl-containing pharmaceutical compounds are sharp, symmetrical and easily integrated.

Catecholamines

CHROM
8903

1. Noradrenaline (norepinephrine)
2. Adrenaline (epinephrine)
3. Dopamine (3-Hydroxytyramine)
4. 3,4-Dihydroxyphenylalanine (DOPA)
5. Phenylalanine
6. Homovanillic Acid (HVA)

Standard Platinum C18 Column




Column: Platinum C18, 3µm, 53 x 7mm Rocket™
Mobile Phase: A: 0.15% TFA, B: Acetonitrile
Gradient: Time: 0 | 3 | 12 | 15
 %B | 3 | 3 | 65 | 65
Flow Rate: 2.5mL/min
Detector: ELSD

Barbiturates

CHROM
8684

Platinum EPS C18 Column



1. Barbitol
2. Butalbital
3. Aprobarbital
4. Talbutal
5. Butobarbital
6. Phenobarbital
7. Amobarbital
8. Mephobarbital

Column: Platinum EPS C18, 3µm, 53 x 7mm Rocket
Mobile Phase: 0.010M Sodium Acetate, pH 4.0: Acetonitrile (75:25)
Flow Rate: 3.0mL/min
Detector: UV at 230nm

The GRACE Davison Product Lines

Alltech

DAVISIL

CHROM

JONES
CHROMATOGRAPHY

MODcol

VYDAC


Reverse Elution Order with Standard Platinum™ and Platinum EPS Columns

Separation improves when minor components elute before, rather than after, closely retained major components. The Standard Platinum / Platinum EPS Column combination lets you choose elution order, which can improve separation — another powerful method development tool.

Vasodilators

CHROM 8387

Standard Platinum C18 Column




NOTE

The elution order changes between Standard Platinum and Platinum EPS Columns.

Platinum EPS C18 Column

CHROM 8388



- | | | |
|-------------------------|----------------------|---|
| 1. Pindolol | Column: | 5µm, 150 x 4.6mm |
| 2. Pentoxifylline | Mobile Phase: | 0.025M KH ₂ PO ₄ , pH 3.0: Methanol (50:50) |
| 3. Isosuprine | Flow Rate: | 1.0mL/min |
| 4. Nifedipine Degradant | Detector: | UV at 220nm |
| 5. Nifedipine | | |


Platinum EPS Columns Have High Polar Compound Capacity

This is especially important for early eluting polar compounds which are often unresolved on conventional reversed phase columns.

Local Anesthetics

CHROM 8159


Standard Platinum C18 Column



Platinum EPS C18 Column

NOTE

Polar compounds interact strongly with Platinum for increased capacity



- | | | |
|----------------|----------------------|--|
| 1. Procaine | Column: | 5µm, 150 x 4.6mm |
| 2. Mepivacaine | Mobile Phase: | A: 0.025M KH ₂ PO ₄ ,
B: pH 3.0: Methanol (50:50) |
| 3. Tetracaine | Gradient: | Time: 0 12
%B: 35 60 |
| 4. Dibucaine | Flow Rate: | 1.0mL/min |
| | Detector: | UV at 220nm |

Capacity Factor Values (k')

Compound	Standard Platinum C18	Platinum EPS C18
1. Procaine	0.44	0.86
2. Mepivacaine	0.70	1.40
3. Tetracaine	3.04	4.30
4. Dibucaine	4.74	6.28

The GRACE Davison Product Lines

Alltech

DAVISIL

CHROM

JONES
CHROMATOGRAPHY

MODcol



VYDAC

Polar Selectivity Benefits

Use Platinum™ EPS Columns to Reduce Run Times

Reduced retention, better peak shape, and greater sensitivity of dioctylphthalate on the Platinum EPS Column compared to the conventional reversed phase C18 column.

Phthalates




1. Dimethylphthalate
2. Dibutylphthalate
3. Dioctylphthalate

Columns: Rocket™ Column, 3µm, 53 x 7mm
Mobile Phase: Acetonitrile:Water (65:35)
Flow Rate: 2.3mL/min
Detector: UV at 254nm

Sharp, Symmetrical Peaks

Platinum Columns produce sharp, symmetrical peaks for these sensitive metal chelators. The competitive column, despite a higher bonded phase loading, cannot make that claim.

Metal Chelators




NOTE

Platinum Column produces sharp symmetrical peaks for chelators.

NOTE

Chelators interact with metals on competitive reversed-phase packings.



1. 2,4-Pentanedione
2. 2,2-Dipyridyl
3. 8-Hydroxyquinoline

Columns: 5µm, 150 x 4.6mm
Mobile Phase: Methanol:Tetrahydrofuran: 0.025M K₂HPO₄, pH 7 (20:20:60)
Flow Rate: 0.7mL/min
Detector: UV at 254nm

The GRACE Davison Product Lines

Alltech

DAVISIL

CHROM

JONES
CHROMATOGRAPHY

MODcol

VYDAC

Reproducibility that Stands the Test of Time

Reproducible columns are the key to rugged, reproducible methods. Platinum™ Columns meet this goal with reproducible silica synthesis, bonding, and column packing.


Platinum Column's manufacturing process controls selectivity, capacity and efficiency — the three factors that control resolution in your method!

Reproducible Silica Synthesis and Bonding Minimize Capacity and Selectivity Variations

Platinum silica's exceptional consistency is the foundation of Platinum column performance. Stringent quality control standards produce unmatched batch-to-batch reproducibility.

Platinum Silica's Surface Area, Pore Volume and Pore Diameter are Exceptionally Consistent

■ Surface Area	Mean = 207	RSD = 1.6%
● Pore Volume	Mean = 0.52	RSD = 0.64%
◆ Pore Diameter	Mean = 99.6	RSD = 1.5%




Reproducible Column Packing Methods Control Efficiency, Symmetry, and Operating Pressure

Alltech's advanced column packing methods minimize column-to-column variation. Efficiency, symmetry, and operating pressure are all tightly controlled.

Platinum HPLC Columns Are Among the Most Reproducible in the Industry

Capacity Factors for Platinum EPS C18 Columns

◆ Biphenyl	■ Toluene	▲ Cinnamic Acid
● Amitriptyline	× Nitrobenzene	



The GRACE Davison Product Lines



DAVISIL



JONES CHROMATOGRAPHY

MODcol

VYDAC



Stability for Long Column Life

Platinum™ Columns deliver long life and stable performance, even under extreme conditions. Methods developed today will work well tomorrow, next week, and next year.

Extreme pH and Temperature Conditions

Platinum Columns tolerate highly acidic conditions (pH 1, 60°C) for more than 500 hours without losing performance.

Performance Remains Stable, Even After 500+ Hours Under Extreme Conditions (pH 1, 60°C)



1. Uracyl
2. Phenol
3. o-Nitroaniline
4. Toluene
5. 4-Hexylbenzoic Acid
6. 4-Octylbenzoic Acid


Column: 5µm, 150 x 4.6mm
Mobile Phase: Acetonitrile:
 0.05M KH₂PO₄, pH 3.2
 (65:35)
Flow Rate: 1.0mL/min
Detector: UV at 254nm

Column Capacity and Selectivity Remain Constant


Platinum Columns do more than just withstand low-pH and high-temperature conditions; they deliver constant capacity and selectivity.

Capacity Factors are Stable Throughout the Study

Standard Platinum C18 Columns



Platinum EPS C18 Columns



The GRACE Davison Product Lines

Alltech

DAVISIL

CHROM

JONES
CHROMATOGRAPHY

MODcol


VYDAC

Platinum™ Columns Simplify Method Development

Use Standard Platinum Columns for Neutral and Moderate Polarity Compounds

Standard Platinum Columns offer less silanol interaction than Platinum EPS Columns. Neutral to moderately polar compounds with few polar functional groups interact weakly with controlled silanols and are retained primarily by reversed-phase mechanisms. For example, nitroanilines are retained and separated on Standard Platinum Columns primarily by hydrophobic interactions.

Positional Isomers of Nitroaniline




Use Platinum EPS Columns for Greater Separation of Compounds with Several Polar Groups

Highly polar compounds with more than two polar functional groups interact strongly with Platinum's silanols. The degree of polar interaction varies by functional group as follows:

NH₂ > COOH > CHO > CN =OH > NO₂ > OCH₃ >
Cl > F > Br > CH₃

Fluphenazine



The GRACE Davison Product Lines

Alltech

DAVISIL

CHROM

JONES
CHROMATOGRAPHY

MODcol

VYDAC


Platinum™ Columns Solve Tough Pharmaceutical Analysis Problems

In pharmaceutical analysis, Platinum Column's utility ranges from general USP methods to specific in-process testing. In drug purity analysis, Platinum's polar selectivity increases resolution between parent compounds and degradants. During method development, Standard Platinum and Platinum EPS Columns can be used to generate alternative selectivities to most reversed-phase columns.

Benzodiazepines and Metabolites

CHROM
8691

- | | |
|-------------------------------|-------------------------------|
| 1. Chlordiazepoxide degradant | 5. Chlordiazepoxide |
| 2. Nitrazepam | 6. Chlordiazepoxide degradant |
| 3. Norchlordiazepoxide | 7. Diazepam |
| 4. Nordiazepam | |




Column: Platinum EPS C18, 3µm, 53 x 7mm Rocket™
Mobile Phase: 0.05M Ammonium Acetate, pH5.5: Acetonitrile (65:35)
Flow Rate: 3.0mL/min
Detector: UV at 254nm

Antibacterials

CHROM
8711

1. Sulfanilimide
2. Cefaclor
3. Cefatrizine
4. Thiamphenicol
5. Cefotaxime




Column: Platinum EPS C18, 5µm, 150 x 4.6mm
Mobile Phase: 0.025M KH₂PO₄, pH 3: Acetonitrile (90:10)
Flow Rate: 1.0mL/min
Detector: UV at 230nm

Famotidine

CHROM
8758

1. Famotidine
2. Sodium Benzoate
3. Methylparaben
4. Propylparaben



Column: Platinum EPS C18, 5µm, 150 x 4.6mm
Mobile Phase: A: 50mM KH₂PO₄, pH 3.0 B: ACN


Time:	0	5	15	20
%B:	10	10	40	40

Flow Rate: 1.0mL/min
Detector: UV at 254nm

Biogenic Amines

CHROM
8297

1. Histamine
2. Tyramine
3. Serotonin
4. β-Phenylethylamine
5. Tryptamine




Column: Platinum EPS C18, 5µm, 250 x 4.6mm
Mobile Phase: Acetonitrile:0.03M KH₂PO₄ pH 3.2 (20:80)
Flow Rate: 1.0mL/min
Detector: UV at 220nm

The GRACE Davison Product Lines

Sedatives/Hypnotics

CHROM
8160

1. Barbitol
2. Phenobarbital
3. Butalbital
4. Oxazepam
5. Triazolam




Column: Standard Platinum™ C18, 5µm, 150 x 4.6mm
Mobile Phase: Methanol:Acetonitrile:0.05M NH₄H₂PO₄, pH 4.5, (50:10:40)
Flow Rate: 1.0mL/min
Detector: UV at 230nm

Central and Respiratory Stimulants

CHROM
8169

1. Theophylline
2. Caffeine
3. Ethamivan
4. 8-Chlorotheophylline




Column: Platinum EPS C18, 5µm, 150 x 4.6mm
Mobile Phase: Methanol:Acetonitrile:1% Acetic Acid, pH 3.2 (30:10:60)
Flow Rate: 1.0mL/min
Detector: UV at 254nm

Anti-Convulsants

CHROM
8163

1. Ethotoin
2. Methsuximide
3. Methetoin
4. Carbamazepine




Column: Platinum EPS C18, 5µm, 150 x 4.6mm
Mobile Phase: Methanol:Acetonitrile:0.005M NaH₂PO₄, pH=5.5, (30:20:50)
Flow Rate: 1.0mL/min
Detector: UV at 220nm

Thyroid Chemicals

CHROM
8282

1. DL-Tyrosine
2. 3-Iodo-Tyrosine
3. 3,5- Diiodo-D-Tyrosine
4. DL-Thyroxine
5. 3,3,5-Triiodo-L-Thyronine



Column: Platinum EPS C18, 5µm, 150 x 4.6mm
Mobile Phase: 0.05M KH₂PO₄, pH 3.2:Acetonitrile (65:35)
Flow Rate: 1.0mL/min
Detector: UV at 254nm

The **GRACE Davison** Product Lines

Alltech

DAVISIL

CHROM

JONES
CHROMATOGRAPHY

MODcol

VYDAC




Extended Polar Selectivity Separates Environmental Contaminants

Platinum™ Column's extended polar selectivity easily resolves closely-related environmental contaminants for improved detection and quantitation. Conventional HPLC columns frequently cannot separate structurally-similar toxic environmental contaminants. Platinum Columns separate triazine herbicides, phenols, and pesticides with exceptional resolution.

Triazine Herbicides

CHROM-8671

1. Simazine
2. Simetryn
3. Prometon
4. Atrazine
5. Ametryn
6. Prometryn
7. Terbutryn




Column: Platinum EPS C18, 1.5µm, 33 x 7mm Rocket™
Mobile Phase: 0.025M KH₂PO₄ pH3:Acetonitrile (65:35)
Flow Rate: 2.0mL/min
Detector: UV at 254nm

Organochlorine Pesticides

CHROM-8437

1. Methoxychlor
2. Dieldrin
3. Heptachlor
4. 4,4-DDT




Column: Standard Platinum C18, 5µm, 150 x 4.6mm
Mobile Phase: Methanol:0.025M K₂HPO₄, pH 7 (75:25)
Flow Rate: 1.0mL/min
Detector: UV at 220nm

Carbamate and Urea Pesticides

CHROM-8293

- | | |
|----------------|-------------|
| 1. Oxamyl | 8. Siduron |
| 2. Methomyl | 9. Diuron |
| 3. Fenuron | 10. Linuron |
| 4. Aldicarb | 11. Neburon |
| 5. Monuron | |
| 6. Fluometuron | |
| 7. Methiocarb | |



Column: Platinum EPS C18, 5µm 150 x 4.6mm
Mobile Phase: **A:** 0.025M KH₂PO₄, pH 3.2
B: Acetonitrile
Gradient:


Time:	0	2	20	30
%B:	20	20	60	60

Flow Rate: 1.0mL/min
Detector: UV at 240nm

Phenols

CHROM-8289

- | | |
|----------------------|----------------------------|
| 1. Phenol | 5. 2-Nitrophenol |
| 2. 4-Nitrophenol | 6. 4-Chloro-3-methylphenol |
| 3. 2,4-Dinitrophenol | 7. 2,4-Dichlorophenol |
| 4. 2-Chlorophenol | 8. Trichlorophenol |
| | 9. Pentachlorophenol |



Column: Standard Platinum C18, 5µm 150 x 4.6mm
Mobile Phase: **A:** 1% Glacial Acetic Acid in Water
B: 1% Glacial Acetic Acid in Methanol
Gradient:

Time:	0	20
%B:	35	95

Flow Rate: 1.0mL/min
Detector: UV at 280nm


The GRACE Davison Product Lines

Platinum™ Column's Dual Retention Mechanism Streamlines Food Analysis

Platinum Columns makes critical food QC methods easier to develop and simpler and faster to run. One Platinum Column, with both reversed-phase and silanophilic retention mechanisms, often replaces two conventional columns in complex food analysis methods.

Catechins in Green Tea


CHROM
8825



Column: Platinum EPS C18, 1.5µm, 33 x 7mm Rocket™
Mobile Phase: 0.05% v/v TFA in 87:13 Water:ACN
Flow Rate: 2.5mL/min
Detector: UV at 210nm

Fat Soluble Vitamins


CHROM
8280



Column: Platinum EPS C18, 5µm, 150 x 4.6mm
Mobile Phase: Acetonitrile:Methanol (95:5)
Flow Rate: 1.0mL/min
Detector: UV at 280nm

Organic Acids


CHROM
8244



Column: Platinum EPS C8, 5µm, 150 x 4.6mm
Mobile Phase: 0.025M KH₂PO₄, pH 2.5:
Methanol (97:3)
Flow Rate: 0.7mL/min
Detector: UV at 220nm

Water Soluble Vitamins

CHROM
8154



Column: Standard Platinum C18, 5µm, 150 x 4.6mm
Mobile Phase: **A:** 0.1M NH₄H₂PO₄ Buffer, pH4.5
B: Acetonitrile
Gradient:

Time:	0	4	15	18
%B:	5	5	40	40

Flow Rate: 1.0mL/min
Detector: UV at 254nm

The **GRACE Davison** Product Lines

Alltech

DAVISIL

CHROM

JONES
CHROMATOGRAPHY

MODcol


VYDAC

Platinum™ EPS Columns Easily Separate Highly Polar Biotech Products

Platinum Column's dual-mode reversed-phase/silanophilic retention mechanism is a perfect match for biotechnology applications. Biopolymers and their building blocks all have structures with polar and non-polar regions. Platinum packing's bonded phase interacts with the non-polar regions, while the controlled silanols interact with the polar regions. This powerful combination is a versatile tool for biotechnology applications.

Tryptic Digest of Bovine Serum Albumin


CHROM
8852



Column: Standard Platinum C18, 300Å, 5µm, 150 x 4.6mm
Mobile Phase: **A:** 0.15% TFA in Water
B: 0.13% TFA in 95:5 Acetonitrile:Water
Gradient: **Time:** 0 | 60 |
%B: 5 | 35 |
Flow Rate: 1.0mL/min
Detector: UV at 216nm

Nucleosides

CHROM
8245




1. Cytidine
2. Uridine
3. Xanthine
4. Guanosine
5. Adenosine

Column: Platinum EPS C18, 5µm, 150 x 4.6mm
Mobile Phase: **A:** 0.03M KH₂PO₄, pH 3.2
B: Acetonitrile
Gradient: **Time:** 0 | 2 | 20 |
%B: 5 | 5 | 30 |
Flow Rate: 0.7mL/min
Detector: UV at 260nm

Peptides (≥5 Amino Acid Residues)

CHROM
8723




1. Peptide T
2. (Arg-8)-Vasotocin
3. (Lys-8)-Vasopressin
4. (Arg-8)-Vasopressin
5. Methionine Enkephalin
6. Oxytocin
7. Leucine Enkephalin
8. Physalaemin
9. Bombesin
10. Angiotensin I

Column: Platinum EPS C18, 5µm, 150 x 4.6mm
Mobile Phase: **A:** 0.15% v/v TFA in Water
B: 0.13% v/v TFA in 95:5 Acetonitrile:Water
Gradient: **Time:** 0 | 2 | 30 |
%B: 15 | 15 | 55 |
Flow Rate: 1.0mL/min
Detector: UV at 216nm

Monophosphate Nucleotides

CHROM
8296



1. Cytidine-5'-Monophosphate
2. Uridine-5'-Monophosphate
3. Guanosine-5'-Monophosphate
4. Adenosine-5'-Monophosphate

Column: Platinum EPS C18, 5µm, 150 x 4.6mm
Mobile Phase: Methanol:0.05M Ammonium Acetate pH 5.5 (5:95)
Flow Rate: 1.0mL/min
Detector: UV at 254nm

The GRACE Davison Product Lines

Platinum™ HPLC Columns

Platinum HPLC Columns are available in a variety of formats. Choose from Analytical and Solvent reducing formats to preparative formats. Custom sizes are also available.



High-Speed Rocket™ Columns Reduce Analysis Time and Solvent Consumption

Separations normally performed on conventional 150mm x 4.6mm columns run faster on Rocket Columns and use less mobile phase, saving both time and money.



Actual Size (33 x 7mm)

Standard Platinum/ Platinum EPS Media Specifications

	100Å Packings
Base Material:	Spherical Silica
Particle Size:	1.5, 3, 5 and 10µm
Pore Volume:	0.51cc/g
Surface Area:	200m ² /g

PLATINUM ALL-GUARD™ CARTRIDGES*

PACKING	I.D. X LENGTH	QTY	PART No.
C18	3.0 x 7.5mm	3	99115
	4.6 x 7.5mm	3	32606
EPS C18	3.0 x 7.5mm	3	99117
	4.6 x 7.5mm	3	32607
C8	4.6 x 7.5mm	3	32612
EPS C8	4.6 x 7.5mm	3	32614
Phenyl	4.6 x 7.5mm	3	32619
Cyano	4.6 x 7.5mm	3	32620
Amino	4.6 x 7.5mm	3	32621
Silica	4.6 x 7.5mm	3	32622
SAX	4.6 x 7.5mm	3	32787
All-Guard Cartridge Holder (Includes Direct-Connect™ Column Coupler)		ea	80101

*Guard Holder required

PLATINUM PREP-GUARD COLUMNS*

PACKING	I.D. X LENGTH	QTY	PART No.
EPS C18	7 x 33mm	1	970216
Direct-Connect Column Coupler		ea	28195

*Direct-Connect Column Coupler required

PLATINUM HPLC COLUMNS			
PACKING	FORMAT	I.D. X LENGTH	PART No.
C18, 1.5µm	Rocket	7 x 33mm	50527
	Rocket	7 x 53mm	50529
C18, 3µm	Solvent Reducer	3.0 x 150mm	32794
	Analytical	4.6 x 100mm	32007
	Analytical	4.6 x 150mm	32020
	Rocket	7 x 33mm	50525
	Rocket	7 x 53mm	50523
C18, 5µm	Solvent Reducer	3.0 x 150mm	32793
	Solvent Reducer	3.0 x 250mm	32792
	Analytical	4.6 x 150mm	32043
	Analytical	4.6 x 250mm	32064
EPS C18, 1.5µm	Rocket	7 x 33mm	50577
	Rocket	7 x 53mm	50579
EPS C18, 3µm	Solvent Reducer	3.0 x 150mm	32799
	Analytical	4.6 x 100mm	32158
	Analytical	4.6 x 150mm	32183
	Rocket	7 x 33mm	50575
	Rocket	7 x 53mm	50573
EPS C18, 5µm	Solvent Reducer	3.0 x 150mm	32806
	Solvent Reducer	3.0 x 250mm	32802
	Analytical	4.6 x 150mm	32214
	Analytical	4.6 x 250mm	32246
	Prep	10 x 150mm	82083
	Prep	22 x 150mm	82093
	Prep	22 x 250mm	82095
C8, 1.5µm	Rocket	7 x 53mm	50529
C8, 3µm	Rocket	7 x 33mm	50532
C8, 5µm	Analytical	4.6 x 150mm	32370
	Analytical	4.6 x 250mm	32375
EPS C8, 3µm	Analytical	4.6 x 250mm	32943
	Analytical	4.6 x 150mm	32415
	Rocket	7 x 33mm	50583
EPS C8, 5µm	Rocket	7 x 53mm	50585
	Analytical	4.6 x 150mm	32420
Phenyl, 3µm	Analytical	4.6 x 250mm	32425
	Analytical	4.6 x 150mm	32631
Phenyl, 5µm	Analytical	4.6 x 150mm	32636
	Analytical	4.6 x 250mm	32641
Cyano, 3µm	Rocket	7 x 33mm	50593
	Rocket	7 x 53mm	50595
Cyano, 5µm	Analytical	4.6 x 150mm	32672
	Analytical	4.6 x 250mm	32681
Amino, 3µm	Analytical	4.6 x 150mm	32706
	Rocket	7 x 53mm	50545
Amino, 5µm	Analytical	4.6 x 150mm	32713
	Analytical	4.6 x 250mm	32722
Silica, 3µm	Analytical	4.6 x 150mm	32535
Silica, 5µm	Analytical	4.6 x 150mm	32542
	Analytical	4.6 x 250mm	32549
SAX, 3µm	Analytical	4.6 x 150mm	32952
SAX, 5µm	Analytical	4.6 x 150mm	32944
	Analytical	4.6 x 250mm	32943

© Copyright 2005 Alltech Associates, Inc. Alltech Associates Inc., is a wholly owned subsidiary of W. R. Grace & Co.-Conn.

Grace Davison is a business unit of W. R. Grace & Co.-Conn., Global Headquarters 7500 Grace Drive, Columbia, MD 21044

The information contained herein is based on our testing and experience and is offered for the user's consideration, investigation and verification. Since operating and use conditions vary and since we do not control such conditions, we must DISCLAIM ANY WARRANTY, EXPRESSED OR IMPLIED, with regard to results to be obtained from the use of this product. Test methods are available on request. Alltech reserves the right to change prices and/or specifications without prior notification. Printed in the USA.

The **GRACE Davison** Product Lines



DAVISIL



Contact your Alltech office or distributor for current or local prices.

© Copyright 2005 Alltech Associates, Inc.

Stock up and Save with Alltech Start-Up Packages

Basic Start-Up Package

Provides the essential products for HPLC system start-up. No frills and no unnecessary cost.



SAVE 20%*

Advanced Start-Up Package

Keeps your HPLC consumables stocked full with products that maximize system and lab efficiency.



SAVE 25%*

BASIC START-UP PACKAGE			
DESCRIPTION	MATERIAL	QTY.	PART No.
Basic Start-Up Package	—	1	97817
Replacement Individual Components			
Short Hex-Head Nuts/Ferrules	PEEK	10	35744
Long Hex-Head Nuts/Ferrules	PEEK	10	35735
Short Finger-Tight™ One Piece*	PEEK	12	32070
1/8" Flange-Free™ Nuts	PEEK	10	37075
1/8" Flange-Free™ Ferrules	Tefzel®	10	20242
1/16" x 0.020" i.d. Tubing	PEEK	10'	35708
1/16" x 0.007" i.d. Tubing	PEEK	10'	35712
1/8" x 0.063" i.d. Tubing	TFE	10'	20063
Plastic Tubing Cutter	—	1	3206
Replacement Blade	—	1	3214
10µm Solvent Inlet Filter	SS	2	28870
20µL Sample Loop	PEEK	1	32195
Open-End Wrench, 1/4" x 5/16"	—	1	1998
Union	PEEK	1	32141

ADVANCED START-UP PACKAGE			
DESCRIPTION	MATERIAL	QTY.	PART No.
Advanced Start-Up Package	—	1	97818
Replacement Individual Components			
Short Hex-Head Nuts/Ferrules	PEEK	10	35744
Long Hex-Head Nuts/Ferrules	PEEK	10	35735
Short Finger-Tight™ No-Slip Fitting*	PEEK	24	32937
Long Finger-Tight™ No-Slip Fitting	PEEK	10	32295
1/8" Flange-Free™ Nuts	PEEK	10	37075
1/8" Flange-Free™ Ferrules	Kel-F®	10	201251
1/16" x 0.020" i.d. Tubing	PEEK	20'	35708
1/16" x 0.005" i.d. Tubing	PEEK	20'	35714
1/8" x 0.063" i.d. Tubing	TFE	10'	20063
Clean-Cut™ Tubing Cutter	—	1	35902
Replacement Blade	—	2	35903
10µm Last-Drop™ Inlet Filter	PAT®	2	28839
In-Line Filter, 2µm	PAT®	1	68250
Replacement Filter Elements	PAT®	5	68152
20µL Flex-Connect™ Sample Loop	PEEK	1	32195
100µL Flex-Connect™ Sample Loop	PEEK	1	32204
Handy-Lok Fitting Wrench	—	1	35751
Open-End Wrench, 1/4" x 5/16"	—	1	1998
Union	PEEK	2	32141

SAVE * Purchasing a package saves you the indicated percentage compared to individual component list prices.

The trademarks listed below and referred to herein are the property of Alltech Associates, Inc.
 ALL-GUARD™ FINGER-TIGHT™ JONES CHROMATOGRAPHY™ ROCKET™
 ALLTECH™ FLANGE-FREE™ LAST-DROP™ VYDAC®
 DAVISIL™ FLEX-CONNECT™ MODCOL®
 DIRECT-CONNECT™ GROM™ PLATINUM™

The trademarks listed below and referred to herein are the property of the corresponding listed entities.
 CLEAN-CUT™ VICI Jour
 KEL-F® Minnesota Mining and Manufacturing Company
 PAT® Systec Inc.
 TEFZEL® E.I. DuPont DeNemours and Company

GRACE®, GRACE DAVISON®, DAVISIL® are registered trademarks of W.R. Grace & Co. - Conn.

This list has been compiled using available published information as of the publication date of this brochure and may not accurately reflect current trademark ownership.

Visit www.chromatography.com for the latest in technical information and applications.

Contact your Alltech office or distributor for current or local prices.

GRACE Davison

2051 Waukegan Rd • Deerfield, IL 60015
 Tel: 847-948-8600 • Fax: 847-948-0477

The GRACE Davison Product Lines

07/05 #444A



DAVISIL



MODCOL



Contact your Alltech office or distributor for current or local prices.

© Copyright 2005 Alltech Associates, Inc.